

## Human NOL3 Recombinant Protein

Catalog Number: PROTO60936

### Overview

Product Name	Human NOL3 Recombinant Protein
Description	Human NOL3 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-208aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 10% glycerol, 1mM DTT, 1mM EDTA

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, O60936, 1-208aa; MGNAQER PSETIDRERK RLVETLQADS GLLLDALLAR GVLTGPEYEA LDALPDAERR VRRLLLLVQG KGEAACQELL  
RCAQRTAGAP DPAWDWQHVG PGYRDRSYDP PCPGHWTP EA PGSGTTCPL PRASDPDEAG GPEGSEAVQS GTPPEPEPEL EAEASKEAEP  
EPEPEPELEP EAEAEPEPEL EPEPDPEPEP DFEERDESED S

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human NOL3 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.