

Human VPS26A Recombinant Protein

Catalog Number: PROTO75436

Overview

Product Name	Human VPS26A Recombinant Protein
Description	Human VPS26A Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-327aa. Application(s): SDS-PAGE.
Size	Starting from 4ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 0.15M NaCl

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, O75436, 1-327aa; MGSSHHHHHH SGLVPRGSH MGSMSFLGGF FGPICEIDIV LNDGETRKMA EMKTEDGKVE KHYLFYDGES
VSGKVNLA FK QPGKRLEHQG IRIEFVGQIE LFNDKSN THE FVNLVKELAL PGELTQSR SY DFEFMQVEKP YESYIGANVR LRYFLKVTIV
RRLTDLVKEY DLIVHQLATY PDVNNSIKME VGIEDCLHIE FEYNKSKYHL KDVIVGKIYF LLVRIKIQHM ELQLIKKEIT GIGPSTTTET ETIAKYEIMD
GAPVKGESIP IRLFLAGYDP TPTMRDVNKK FSVRYFLNLV LVDEEDRRYF KQQEII LRWK APEKLRKQRT NFHQRFESPE SQASAEQPEM

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human VPS26A Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.