

Mouse Agr2 Recombinant Protein

Catalog Number: PROTO88312

Overview

Product Name	Mouse Agr2 Recombinant Protein
Description	Mouse Agr2 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 21-175aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 20% glycerol, 1mM DTT, 1mM EDTA

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Mouse, O88312, 21-175aa; MGSSHHHHHH SGLVPRGSH MGSHMKDTTV KSGAKKDPKD SRPKLPQTLS RGWGDQLIWT QTYEEALYRS
KTSNRPLMVI HHLDECPHSQ ALKKVFAEHK EIQLAEQFV LLNLVYETTD KHLSPDGQYV PRIVFVDP SL TVRADITGRY SNRLYAYEPS
DTALLYDNMK KALKLLKTEL

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Agr2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.