

## Mouse Legumain/Asparaginyl Endopeptidase Recombinant Protein

Catalog Number: PROTO89017

### Overview

Product Name	Mouse Legumain/Asparaginyl Endopeptidase Recombinant Protein
Description	Mouse Legumain/Asparaginyl Endopeptidase Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 18-435aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

0.5mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Mouse, O89017, 18-435aa; VPGVDDDPED GGKHWVVIVA GSNGWYNYRH QADACHAYQI IHRNGIPDEQ IIVMMYDDIA NSEENPTPGV  
VINRPNGTDV YKGV LKDYTG EDVTPENFLA VLRGDAEAVK GKGSGKVLKS GPRDHVFIYF TDHGATGILV FPNDLHVVD LNKTI RYMYE  
HKMYQKMFY IEACESGSMM NHLPDDIN VY ATTAANPKES SYACYDEER GTYLGDWYSV NWMEDSDVED LTKETLHKQY HLVKSHTNTS  
HVMQYGNKSI STMKVMQFQG MKHRASSPIS LPPVTHLDTL PSPDVPLTIL KRKLLRTNDV KESQNLIGQI QQFLDARHVI EKSVHKIVSL  
LAGFGETAER HLSERTMLTA HDCYQEA VTH FRTHCFNWHS VTYEHALR YL V LANC EAP YPIDRIEMAM DKVCLSHYLE HHHHHH

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Legumain/Asparaginyl Endopeptidase Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.