

CLSTN1 (NM_001009566) Human Recombinant Protein

Catalog Number: PROTO94985

Overview

Product Name	CLSTN1 (NM_001009566) Human Recombinant Protein
Description	Recombinant protein of human calyntenin 1 (CLSTN1), transcript variant 1
Size	20 µg
Tag	C-Myc/DDK
Form	Frozen Solution in PBS Buffer
Source	HEK293T
Formulation	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

Concentration

>50 ug/mL as determined by microplate BCA method

Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Amino Acid Sequence

MLRRPAPALAPAARLLL LAGLLCGGGVWAARV NKHKPWLEPTYHGIVTENDNTVLLDPPLIALDKDAPLRF AEFVTVTKEGEICGFKIHGQNV PFDVAVVDKSTGEGVIRSKEKLDCE LQKDY SFTIQAYDCGKGP DGTNVK KSHKATVHIQVNDVNEYAPVFKEKSYKATVIEGKQYDSILRVEAVDADCSPQFSQICSYEIITPDVPFTVDKDG YIKNTEKLN YGKEHQYKLT V TAYDCGKKRATEDVLVKISIKPTCTPGWQGWNNRIEYEPGTGALAVFPNIHLETCD E P V A S V Q A T V E L E T S H I G K G C D R D T Y S E K S L H R L C G A A A G T A E L L P S P S G S L N W T M G L P T D N G H D S D Q V F E F N G T Q A V R I P D G V V S V S P K E P F T I S V W M R H G P F G R K K E T I L C S S D K T D M N R H H Y S L Y V H G C R L I F L R Q D P S E E K K Y R P A E F H W K L N Q V C D E E W H H Y V L N V E F P S V T L Y V D G T S H E P F S V T E D Y P L H P S K I E T Q L V V G A C W Q E F S G V E N D N E T E P V T V A S A G G D L H M T Q F F R G N L A G L T L R S G K L A D K K V I D C L Y T C K E G L D L Q V L E D S G R G V Q I Q A H P S Q L V L T L E G E D L G E L D K A M Q H I S Y L N S R Q F P T P G I R R L K I T S T I K C F N E A T C I S V P P V D G Y V M V L Q P E E P K I S L S G V H H F A R A A S E F E S E G V F L F P E L R I I S T I T R E V E P E G D G A E D P T V Q E S L V S E E I V H D L D T C E V T V E G E E L N H E Q E S L E V D M A R L Q Q K G I E V S S E L G M T F T G V D T M A S Y E E V L H L L R Y R N W H A R S L L D R K F K L I C S E L N G R Y I S N E F K V E V N V I H T A N P M E H A N H M A A Q P Q F V H P E H R S F V D L S G H N L A N P H P A V V P S T A T V V I V V C V S F L V F M I I L G V F R I R A A H R R T M R D Q D T G K E N E M D W D D S A L T I T V N P M E T Y E D Q H S S E E E E E E E E E E S E D G E E E D D I T S A E S E S E E E E G E Q G D P Q N A T R Q Q Q L E W D D S T L S Y

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



CLSTN1 (NM_001009566) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.