

Human TL1A/TNFSF15 Recombinant Protein

Catalog Number: PROTO95150

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human TL1A/TNFSF15 Recombinant Protein
Description	Human TL1A/TNFSF15 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 72-251aa. Application(s): SDS-PAGE, Bioactivity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.5) containing 0.15M NaCl, 40% glycerol

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, O95150, 72-251aa; LKGQEFAPS HQQVYAPLRA DGDKPRAHLT VVRQTPTQHF KNQFPALHWE HELGLAFTKN RMNYTNKFLI
IPESGDYFIY SQVTRGMTS ECSEIRQAGR PNKPSITVV ITKVTDSTPE PTQLLMGTS VCEVGSNWVQ PIYLGAMFSL QEGDKLMVNV
SDISLVDYTK EDKTFFGAFL L

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human TL1A/TNFSF15 Recombinant Protein](#)

For Research Use Only. Not for use in diagnostic procedures.