

Human STBD1 Recombinant Protein

Catalog Number: PROTO95210

Overview

Product Name	Human STBD1 Recombinant Protein
Description	Human STBD1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 24-358aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% glycerol, 1mM DTT

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, O95210, 24-358aa; MGSSHHHHHH SGLVPRGSH MGSRRGGPDT GKDGDAEQEK DAPLGAAIP GGHQSGSSGL SPGPSSGQELV
TKPEHLQESN GHLISKTKDL GKLQAASWRL QNPSREVCDN SREHVPSGQF PDTEAPATSE TSNSRSYSEV SRNESLESPM GEWGFQKGQE
ISAKAATCFA EKLPSNLLK NRAKEEMSL S DLNSQDRVDH EEWEMVPRHS SWGDVGVGGS LKAPVLNLNQ GMDNGRSTLV EARGQQVHGK
MERVAVMPAG SQQVSVRFQV HYVTSTDVQF IAVTGDHECL GRWNTYIPLH YNKDGFWSHS IFLPADTVVE WKFVLVENG VTRWEECSNR
FLETGHEDKV VHAWWGIH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human STBD1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.