

Human DcR3/TNFRSF6B Recombinant Protein

Catalog Number: PROTO95407

Overview

Product Name	Human DcR3/TNFRSF6B Recombinant Protein
Description	Human DcR3/TNFRSF6B Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 30-300aa. Application(s): SDS-PAGE, Denatured. Denatured protein is most likely not the best option for functional studies. It is better suited for Western Blot (WB) or imaging assays.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol 0.4M urea

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, O95407, 30-300aa; VAETPTY PWRDAETGER LVCAQCPPGT FVQRPCRDS PTTGCPGPPR HYTQFWNYLE RCRYCNVLCG
EREEARACH ATHNRACRCR TGFFAHAGFC LEHASCPPGA GVIAPGTPSQ NTQCQPCPPG TFSASSSSSE QCQPHRNCTA LGLALNVPGS
SSHDTLCTSC TGFPLSTRVP GAEECERAVI DFVAFQDISI KRLQRLQLAL EAPEGWGPTP RAGRAALQLK LRRRLTELLG AQDGALLVRL
LQALRVARMP GLERSVRERF LPVH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human DcR3/TNFRSF6B Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.