

Human FGFR1OP Recombinant Protein

Catalog Number: PROTO95684

Overview

Product Name	Human FGFR1OP Recombinant Protein
Description	Human FGFR1OP Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-379aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 1mM DTT

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, O95684, 1-379aa; MGSSHHHHHH SGLVPRGSH MGSMAATAA AVVAEEDTEL RDLLVQTLEN SGVLNRIKAE LRAAVFLALE
EQEKVENKTP LVNESLKKFL NTKDGRVLAS LVAEFLQFFN LDFTLAVFQP ETSTLQGLEG RENLARDLGI IEAEGTVGGP LLELVIRRCQ
QKEKGPTTGE GALDLSDVHS PPKSPEGKTS AQTPSKKAN DEANQSDTSV SLSEPKSKSS LHLLSHETKI GSFLSNRTLK GKDKAGLCPD
EDDMEGDSFF DDPIPKPEKT YGLRKEPRKQ AGSLASLSDA PPLKSGLSSL AGAPSLKDSE SKRGNTVLKD LKLISDKIGS LGLGTGEDDD
YVDDFNSTSH RSEKSEISIG EEIEEDLSVE IDDINTSDKL DDLTQDLTVS QLSDVADYLE DVA

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human FGFR1OP Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.