

## Human BAG3 Recombinant Protein

Catalog Number: PROTO95817

### Overview

Product Name	Human BAG3 Recombinant Protein
Description	Human BAG3 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-575aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 1mM EDTA, 0.1mM PMSF, 10% glycerol

### Concentration

0.5mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, O95817, 1-575aa; MSAATHSPMM QVASNGDRD PLPPGWEIKI DPQTGWPFV DHNSRTTWN DPRVPSGPK ETPSSANGPS  
REGSRLPPAR EGHVPVYQLR PGYIPIVLH EGAENRQVHP FHVYPQPGMQ RFRTEAAAAA PQRSQSPLRG MPETTQPKQ CGQVAAAAAA  
QPPASHGPER SQSPAASDCS SSSSSASLPS SGRSSLGSHQ LPRGYISIPV IHEQNVTRPA AQPSFHQAQK THYPAQQGEY QTHQPVYHKI  
QGDDWEPRPL RAASPFRSSV QGASSREGSP ARSSTPLHSP SPIRVHTVVD RPQQPMTHRE TAPVSQPENK PESKPGVGP ELPPGHIPIQ  
VIRKEVDSKP VSQKPPPPSE KVEVKVPPAP VPCPPSPGP SAVPSSPKSV ATEERAAPST APAEATPPKP GEAEAPPKHP GVLKVEAILE  
KVQGLEQAVD NFEGKKTDDK YLMIEEYLT ELLALDSVDP EGRADVRQAR RDGVRKVQTI LEKLEQKAID VPGQVQVYEL QPSNLEADQP  
LQAIMEMGAV AADKGKKNAG NAEDPHTETQ QPEATAAATS NPSSMTDTPG NPAAP

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human BAG3 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.