

Human LILRA3/CD85e Recombinant Protein

Catalog Number: PROTOQ8N6C8

Overview

Product Name	Human LILRA3/CD85e Recombinant Protein
Description	Human LILRA3/CD85e Recombinant Protein expressed in HEK293 with hIgG-His-Tag. Sequence domain: 24-439aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method)
Size	starting from 10ug
Tag	hIgG-His-Tag
Form	Liquid
Source	HEK293
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by Absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino acid Sequence

GPLPKPT LWAEPGSVIT QGSPVTLRCQ GSLETQEYHL YREKKTALWI TRIPQELVKK GQFPILSITW EHAGRYCCIIYGSHTAGLSES SDPLELVVTG
AYSKPTLSAL PSPVVTSGGN VTIQCDSQVA FDGFILCKEG EDEHPQCLNS HSHARGSSRAIFSVGPVSPS RRWSYRCYGY DSRAPYVWSL
PSDLLGLLVP GVSKKPLSV QPGPVVAPGE KLTFCQGS DA GYDRFVLYKEWGRDFLQRPQ RQPQAGLSQA NFTLGPVSRYS YGGQYTCSGA
YNLSSEWSAP SDPLDILITG QIRARPFLSV RPGPTVASGENVTLLCQSQG GMHTFLLTKE GAADSPLRLK SKRQSHKYQA EFPMSPV TSA
HAGTYRCYGS LSSNPYLLTH PSDPLELVVS GAAETLSP PQ NKSDSKAGE

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human LILRA3/CD85e Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.