

Mouse DHFR Recombinant Protein

Catalog Number: PROTP00375

Overview

Product Name	Mouse DHFR Recombinant Protein
Description	Mouse DHFR Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-187aa. Application(s): Enzyme Activity, SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 2mM DTT, 0.1M NaCl.

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Mouse, P00375, 1-187aa; MVRPLNCIVA VSQNMGIGKN GDLWPPLRN EFKYFQRMTT TSSVEGKQNL VIMGRKTWFS IPEKNRPLKD
RINIVLSREL KEPPRGAHFL AKSLDDALRL IEQPELASKV DMVWIVGGSS VYQEAMNQPG HLRLFVTRIM QEFESDTFFP EIDLGKYKLL
PEYPGVLSEV QEEKGIKYKF EVYEKKD

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse DHFR Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.