

e.coli sodA Recombinant Protein

Catalog Number: PROTP00448

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	e.coli sodA Recombinant Protein
Description	E. coli sodA Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-206aa. Application(s): Enzyme Activity, SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% glycerol, 0.1M NaCl

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

E.coli, P00448, 1-206aa; MGSSHHHHHH SGLVPRGSH MSYTLPSLPY AYDALEPHFD KQTMEIHHTK HHQTYVNNAN AALES LPEFA NLPVEELITK LDQLPADKKT VLRNNAGGHA NHSLFWKGLK KGTT LQGD LK AAIERDFGSV DNFKA EF EKA AASRFGSGWA WLVLKGD KLA VVSTANQDSP LMGEAISGAS GFPIMGLDVW EHAYYLKFQN RRPDIKEFW NVVNWDEAAA RFAAKK

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[e.coli sodA Recombinant Protein](#)

For Research Use Only. Not for use in diagnostic procedures.