

Human Ornithine Carbamoyltransferase/OTC Recombinant Protein

Catalog Number: PROTP00480

Overview

Product Name	Human Ornithine Carbamoyltransferase/OTC Recombinant Protein
Description	Human Ornithine Carbamoyltransferase/OTC Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 33-354aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM MES buffer (pH 6.0) containing 2mM DTT, 10% glycerol, 100mM NaCl

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P00480, 33-354aa; NKVQL KGRDLLTKN FTGEEIKYML WLSADLKFRI KQKGEYLP LL QGKSLGMIFE KRSTRTRLST ETGFALLGGH PCFLTQDIH LGVNESLTD T ARVLSSMADA VLARVYKQSD LDTLAKEASI PIINGLS DLY HPIQILADYL TLQEHYSLK GLTLSWIGDG NNILHSIMMS AAKFGMHLQA ATPKGYEPDA SVTKLAEQYA KENGT KLLLT NDPLEAAHGG NVLITDTWIS MGQEEK KKR LQAFQGYQVT MKTAKVAASD WFTLHCLPRK PEEVDDEVFY SPRSLVFPEA ENRKWTIMAV MVSLLDYSP QLQKPKF

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Ornithine Carbamoyltransferase/OTC Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.