

## Human GOT2 Recombinant Protein

Catalog Number: PROTP00505

### Overview

Product Name	Human GOT2 Recombinant Protein
Description	Human GOT2 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 30-430aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% glycerol, 1mM DTT

### Concentration

0.5mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

Human, P00505, 30-430aa; MGSSHHHHHH SGLVPRGSH MGSSSWWTHV EMGPPDPILG VTEAFKRDTN SKKMNLGVGA YRDDNGKPYV  
LPSVRKAEAQ IAAKNLDKEY LPIGGLAEFC KASAEALGE NSEVLKSGRF VTVQTISGTG ALRIGASFLQ RFFKFSRDVF LPKPTWGNHT  
PIFRDAGMQL QGYRYDPKT CGFDFTGAVE DISKIPEQSV LLLHACAHPN TGVDPRPEQW KEIATVVKKR NLFAFFDMAY QGFASGDGDK  
DAWAVRHIE QGINVCLCQS YAKNMGLYGE RVGAFTMVCK DADEAKRVES QLKILIRPMY SNPPLNGARI AAAILNTPDL RKQWLQEVKV  
MADRIIGMRT QLVSNLKKEG STHNWQHITD QIGMFCFTGL KPEQVERLIK EFSIYMTKDG RISVAGVTSS NVGYLAHAIH QVTK

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human GOT2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.