

## Human u-Plasminogen Activator/Urokinase Recombinant Protein

Catalog Number: PROTP00749

### Overview

Product Name	Human u-Plasminogen Activator/Urokinase Recombinant Protein
Description	Human u-Plasminogen Activator/Urokinase Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 21-431aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

0.5mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

Human, P00749, 21-431aa; SNELHQVPSN CDCLNGGTCV SNKYFSNIHW CNCPKKFGGQ HCEIDKSKTC YEGNGHFYRG KASTDTMGRP  
CLPWNSATVL QQTYHAHRSD ALQLGLGKHN YCRNPDNRRR PWCYVQVGLK PLVQECMVHD CADGKKPSSP PEELKFQCGQ KTLRPRFKII  
GGEFTTIENQ PWFAAIYRRH RGGSVTYVCG GSLISPCWVI SATHCFIDYP KKEDYIVYLG RSRLNSNTQG EMKFEVENLI LHKDYSADTL  
AHHNDIALLK IRSKEGRCAQ PSRTIQTICL PSMYNDPQFG TSCEITGFGK ENSTDYLYPE QLKMTVVKLI SHRECQPHY YGSEVTTKML  
CAADPQWKTD SCQGDSSGGL VCSLQGRMTL TGIVSWGRCG ALKDKPGVYT RVSHFLPWIR SHTKEENGLA L

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human u-Plasminogen Activator/Urokinase Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.