

Human recombinant Renin protein

Catalog Number: PROTP00797-5

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human recombinant Renin protein
Description	Human recombinant Renin protein expressed in HEK293 with His-Tag. Sequence domain: 24-406aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	HEK293
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

1mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P00797, 24-406aa; LPTDTTTFKR IFLKRMPSIR ESLKERVDM ARLGPEWSQP MKRLTLGNTT SSVILTNYMD TQYYGEIGIG TPPQTFKVV
DTGSSNVWVP SSKCSRLYTA CVYHKLFDAS DSSSYKHNGT ELTRYSTGT VSGFLSQDII TVGGITVTQM FGEVTEMPAL PFMLAEFDGV
VGMGFIEQAI GRVTPIFDNI ISQGVLKEDV FSFYNNRDE NSQSLGGQIV LGGSDPQHYE GNHFYINLIK TGVWQIQMKG VSVGSSTLLC
EDGCLALVDT GASYISGSTS SIEKLMEALG AKKRLFDYVV KCNEGPTLPD ISFHLGGKEY TLTSADYVFQ ESYSSKKLCT LAIHAMDIPP
PTGPTWALGA TFIRKFYTEF DRRNNRIGFA LAR

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Human recombinant Renin protein

For Research Use Only. Not for use in diagnostic procedures.