

## Human recombinant EGF (Epidermal growth factor) protein, AF

Catalog Number: PROTP01133-5

Product Name	Human recombinant EGF (Epidermal growth factor) protein, AF
Description	Human epidermal Growth Factors (EGF) is a 6 kDa cytokine with 53 amino acid residues. EGF is mainly secreted from ectodermal cells, monocytes, kidney and duodenal glands. Upon binding to its receptor, EGFR, EGF acts to stimulate cell growth and proliferation of epithelial cells, play important roles in many developmental processes including accelerate tooth eruption, inhibits gastric acid secretion, and involve in wound healing.
Size	100ug,500ug,1mg,5mg
Tag	His Tag (C-term)
Form	Lyophilized
Source	Escherichia coli
Formulation	The protein was lyophilized from a 0.2 um filtered solution containing 1X PBS, pH 8.0. If you have any concerns or special requirements, please confirm with us.

### Storage

Lyophilized protein should be stored at -20°C for 1 year. Upon reconstitution, store at 2°C to 8°C for up to 1 week. Further dilute in a buffer containing a carrier protein or stabilizer (e.g. 0.1% BSA, 10%FBS, 5%HSA or 5% trehalose solution), protein aliquots should be stored at -20°C or -80°C for 3-6 months. Avoid repeated freeze/thaw cycles.

### Purity

>95% as determined by SDS-PAGE.

### Molecular weight

The protein has a calculated MW of 7.16 kDa. The protein migrates as 9-11 kDa under reducing condition (SDS-PAGE analysis).

### Activity

Measure by its ability to induce 3T3 cells proliferation. The ED<sub>50</sub> for this effect is 1.0 x 10<sup>2</sup> IU/mg.

### Endotoxin

<0.1 EU per 1 ug of the protein by the LAL method.

### Amino Acid Sequence

MNSDSECLSHDGYCLHDGVCMYIEALDKYACNCVVG YIGERCQYRDLKWWELR with polyhistidine tag at the C-terminus

### Reconstitution

Centrifuge at 3000 rpm for 5 mins before opening. It is recommended to reconstitute the lyophilized protein in sterile H<sub>2</sub>O to a concentration not less than 100 ug/mL and incubate the stock solution at room temperature for at least 20 mins to ensure sufficient re-dissolved. Do Not Vortex! Vigorous shaking may impair the biological activity of the protein.

## Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human recombinant EGF (Epidermal growth factor) protein, AF

For Research Use Only. Not for use in diagnostic procedures.