

Human TSH beta Recombinant Protein

Catalog Number: PROTP01222

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human TSH beta Recombinant Protein
Description	Human TSH beta Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 21-138aa. Application(s): SDS-PAGE, Denatured. Denatured protein is most likely not the best option for functional studies. It is better suited for Western Blot (WB) or imaging assays.
Size	Starting from 4ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.5) containing 10% glycerol, 0.4M Urea

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, P01222, 21-138aa; FCIPT EY TMHIERRECA YCLTINTTIC AGYCMTRDIN GKLF LPKYAL SQDVCTYRDF IYRTVEIPGC PLHVAPYFSY PVALSCKCGK CNTDYSDCIH EAIK TNYCTK PQKSYLVGFS V

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Human TSH beta Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.