

Human recombinant IFN-alpha 1/IFNA1 protein

Catalog Number: PROTP01562-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human recombinant IFN-alpha 1/IFNA1 protein
Description	Human recombinant IFN-alpha 1/IFNA1 protein expressed in E.coli with Non-Tagged. Sequence domain: 24-189aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 20ug
Tag	Non-Tagged
Form	Liquid
Source	E.coli
Formulation	Phosphate-Buffered Saline (pH 7.4)

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P01562, 24-189aa; MCDLPETHSL DNRRTLMLLA QMSRISPSSC LMDRHDFGFP QEEFDGNQFQ KAPAI SVLHE LIQQIFNLFT
TKDSSAAWDE DLLDKFCTEL YQQLNDLEAC VMQEERVGET PLMNADSILA VKKYFRITL YLTEKKYSPC AWEVVRAEIM RSLSLSTNLQ ERLRRKE

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human recombinant IFN-alpha 1/IFNA1 protein](#)

For Research Use Only. Not for use in diagnostic procedures.