

Human Interferon alpha-14 Recombinant Protein

Catalog Number: PROTP01570

Overview

Product Name	Human Interferon alpha-14 Recombinant Protein
Description	Human Interferon alpha-14 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 24-189aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 50% glycerol, 0.2M NaCl, 2mM DTT

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P01570, 24-189aa; MGSSHHHHHH SGLVPRGSH MGSHMCNLSQ THSLNRRRTL MLMAQMRRIS PFSCLKDRHD FEFPQEFDG
NQFQKAQAIS VLHEMMQQTFLFSTKNSSA AWDETLEKF YIELFQQMND LEACVIQEVG VEETPLMNED SILAVKKYFQ RITLYLMEKK
YSPCAWEVVR AEIMRSLSFS TNLQKRLRRK D

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Interferon alpha-14 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.