

Mouse IFN-gamma Recombinant Protein

Catalog Number: PROTP01580

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Mouse IFN-gamma Recombinant Protein
Description	Mouse IFN-gamma Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 23-155aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM MES buffer (pH 5.0) containing 40% glycerol, 0.2M NaCl, 1mM DTT

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Mouse, P01580, 23-155aa; HGTVI ESLESLNNYF NSSGIDVEEK SLFLDIWRNW QKDGMKILQ SQIISFYRLR FEVLKDNQAI SNNISVIESH
LITTFFSNSK AKKDAFMSIA KFEVNNPQVQ RQAFNELIRV VHQLPESSL RKRKRSRC

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Mouse IFN-gamma Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.