

## Mouse IL-3 Recombinant Protein

Catalog Number: PROTP01586

### Overview

|              |  |
|--------------|--|
| Product Name | Mouse IL-3 Recombinant Protein   |
| Description  | Mouse IL-3 Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 27-166aa. Application(s): SDS-PAGE, Bioactivity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method). |
| Size         | Starting from 10ug   |
| Tag          | His-Tag  |
| Form         | Liquid   |
| Source       | Baculovirus  |
| Formulation  | Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol   |

### Concentration

0.25mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Mouse, P01586, 27-166aa; ASISGRDTHR LTRTLNCSSI VKEIIGKLPE PELKTDDEGP SLRNKSFRRV NLSKFVESQG EVDPEDRYVI KSNLQKLNCC LPTSANDSAL PGVFIRLDD FRKKLRFYMV HLNDLETVLT SRPPQPASGS VSPNRTGVEC

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Mouse IL-3 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.