

## Human CD25/IL-2 R alpha Recombinant Protein

Catalog Number: PROTP01589

### Overview

Product Name	Human CD25/IL-2 R alpha Recombinant Protein
Description	Human CD25/IL-2 R alpha Recombinant Protein expressed in Baculovirus with hIgG-His-tag. Sequence domain: 22-240aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hIgG-His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

0.5mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, P01589, 22-240aa; ADPELCDDDDP PEIPHATFKA MAYKEGTM LN CECKRGFRRI KSGSLYMLCT GNSSHSSWDN QCQCTSSATR  
NTTKQVTPQP EEQKERKTE MQSPMQPVDQ ASLPGHCREP PPWENEATER IYHFVVGQMV YYQCVQGYRA LHRGPAESVC KMTHGKTRWT  
QPQLICTGEM ETSQFPGEEK PQASPEGRPE SETSCLVTTT DFQIQTEMAA TMETSIFTTE YQLEPKSCDK THTCPAP PAPER ELLGGPSVFL  
FPPKPKDTLM ISRTPEVTCV VVDVSHEDPE VKFNWYVDGV EVHNAKTKPR EEQYNSTYRV VSVLTVLHQD WLNGKEYKCK VSNKALPAPI  
EKTISKAKGQ PREPQVYTLPSRDELTKNQ VSLTCLVKGF YPSDIAVEWE SNGQPENNYK TTPVLDSDG SFFLYSKLTV DKSRWQQGNV  
FSCVMHEAL HNHYTQKSLSLSPGKHHHHH H

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human CD25/IL-2 R alpha Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.