

## Mouse IgG2A Recombinant Protein

Catalog Number: PROTP01863

### Overview

Product Name	Mouse IgG2A Recombinant Protein
Description	Mouse IgG2A Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 98-330aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Mouse, P01863, 98-330aa; EPRGPTI KPCPPCKCPA PNLLGGPSVF IFPPKIKDVL MISLSPIVTC VVVDVSEDDP DVQISWVFN VEVHTAQTQT HREDYNSTLR VVSALPIQH QDWMSGKEFKC KVNNKDL PAP IERTISKPKG SVRAPQVYVL PPPEEEMTKK QVTLTCMVTD FMPEDIYVEW TNNKTELNY KNTEPVLDS GSYFMYSKLR VEKKNWVERN SYSCSVVHEG LHNHHTTKSF SRTPGK

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Mouse IgG2A Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.