

Human AlphaB Crystallin/CRYAB Recombinant Protein

Catalog Number: PROTP02511

Overview

Product Name	Human AlphaB Crystallin/CRYAB Recombinant Protein
Description	Human AlphaB Crystallin/CRYAB Recombinant Protein expressed in E. coli (non-tagged). Sequence domain: 1-175aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	Non-Tagged
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 7.5) containing 50mM NaCl, 1mM EDTA

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P02511, 1-175aa; MDIAIHHPWI RRPFFPFHSP SRLFDQFFGE HLLSDLFPT STLSPFYLR PPSFLRAPSW FDTGLSEMRL EKDRFSVNLD VKHFSPEELK VKVLGDVIEV HGKHEERQDE HGFISREFHR KYRIPADVDP LTITSSLSSD GVLTVNGPRK QVSGPERTIP ITREEKPAVT AAPKK

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human AlphaB Crystallin/CRYAB Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.