

## Human KRT14 Recombinant Protein

Catalog Number: PROTP02533

### Overview

Product Name	Human KRT14 Recombinant Protein
Description	Human KRT14 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-472aa. Application(s): SDS-PAGE, Denatured. Denatured protein is most likely not the best option for functional studies. It is better suited for Western Blot (WB) or imaging assays.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 0.4M urea

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 85% by SDS-PAGE

### Amino Acid Sequence

Human, P02533, 1-472aa; MTTCSRQFTS SSSMKGSCGI GGGIGGGSSR ISSVLAGGSC RAPSTYGGGL SVSSSRFSSG GAYGLGGGYG GGFSSSSSF  
GSGFGGGYGG GLGTGLGGGF GGGFAGGDGL LVGSEKVTMQ NLNDRLASYL DKVRALEEAN ADLEVKIRDW YQRQRP AEIK DYSPYFKTIE  
DLRNKILTAT VDNANVLLQI DNARLAADD F RTKYETELNL RMSVEADING LRRVLDELTL ARADLEMQIE SLKEELAYLK KNHEEEMNAL  
RGQVGGDVNV EMDAAPGVDL SRILNEMRDQ YEKMAEKNRK DAEWFFTKT EELNREVATN SELVQSGKSE ISELRRTMQN LEIELQSQLS  
MKASLENSLE ETKGRYCMQL AQIQEMIGSV EEQLAQLRCE MEQQNQEYKI LLDVKTRLEQ EIATYRRLL EGEDAHLSSSQ FSSGSQSSRD  
VTSSSRQIRT KVMDVHDGKV VSTHEQVLRT KN

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human KRT14 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.