

## BSA Bovine Serum Albumin protein

Catalog Number: PROTP02769

### Overview

Product Name	BSA Bovine Serum Albumin protein
Description	Contains 583 amino acid residues derived from the prototypical bovine serum albumin sequence. Suitable for use as a biochemical, excipient (an inert substance used as a diluent or vehicle for a protein), in culture media and chromatographic applications.
Size	10gr, 50gr, 100gr
Tag	
Form	Lyophilized freeze dry yellowish powder.
Source	Bovine serum
Formulation	

### Storage

Lyophilized BSA although stable at room temperature for 3 weeks, should be stored 2-8°C. Upon reconstitution BSA should be stored at 4°C between 2-7 days and for future use below -18°C. Please prevent freeze-thaw cycles.

### Purity

Greater than 97.0%.

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



BSA Bovine Serum Albumin protein

For Research Use Only. Not for use in diagnostic procedures.