

## Human AFP Recombinant Protein

Catalog Number: PROTP02771

### Overview

Product Name	Human AFP Recombinant Protein
Description	Human AFP Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 19-609aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

0.5mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

Human, P02771, 19-609aa; RTLHRNE YGIASILDSY QCTAEISLAD LATIFFAQFV QEATYKEVSK MVKDALTAIE KPTGDEQSSG CLENQLPAFL EELCHEKEIL EKYGHSDCCS QSEEGRHNCF LAHKKPTPAS IPLFQVPEPV TSCEAYEEDR ETFMKNKIYE IARRHPFLYA PTILLWAARY DKIIIPSCCKA ENAVECFQTK AATVTKELRE SLLNQHACA VMKNFGTRTF QAITVTKLSQ KFTKVNFTI QKLVLDDVAHV HEHCCRGDVL DCLQDGEKIM SYICSQQDTL SNKITECCKL TTLERGQCII HAENDEKPEG LSPNLNRFLG DRDFNQFSSG EKNIFLASFV HEYSRRHPQL AVSVILRVAK GYQELLEKCF QTENPLECQD KGEELQKYI QESQALAKRS CGLFQKLGEY YLQNAFLVAY TTKAPQLTSS ELMAITRKMA ATAATCCQLS EDKLLACGEG AADIIIGHLC IRHEMTPVNP GVGQCCTSSY ANRRPCFSSL VVDETYVPPA FSDDKFIFHK DLCQAQGVAL QTMKQEFLIN LVKQKPQITE EQLEAVIADF SGLLEKCCQG QEQEVCF AEE GQKLIKTRA ALGV

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human AFP Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.