

## Human FTL Recombinant Protein

Catalog Number: PROTP02792

### Overview

Product Name	Human FTL Recombinant Protein
Description	Human FTL Recombinant Protein expressed in E. coli (non-tagged). Sequence domain: 1-175aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	Non-Tagged
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 7.5)

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, P02792, 1-175aa; MSSQIRQNYSTDVEAAVNSLVNLYLQASYTYLSLGFYFDRDDVALEGVSHFFRELAEEKREGYERLLKMQ  
NQRGGRALFQDIKKPAEDEVGKTPDAMKAA MALEKKNLQA LLDLHALGSA RTDPHLCDFLETHFLDEEVK LIKKMGDHLT NLHRLGGPEA  
GLGEYLFERL TLKHD

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human FTL Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.