

## Human Annexin A1 Recombinant Protein

Catalog Number: PROTP04083

### Overview

Product Name	Human Annexin A1 Recombinant Protein
Description	Human Annexin A1 Recombinant Protein expressed in E. coli (non-tagged). Sequence domain: 1-346aa. Application(s): SDS-PAGE.
Size	Starting from 4ug
Tag	Non-Tagged
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 100mM NaCl, 1mM DTT, 10% glycerol,

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, P04083, 1-346aa; MAMVSEFLKQ AWFIEENEEQE YVQTVKSSKG GPGSAVSPYP TFPSSDVAA LHKAIMVKGV DEATIIDILT  
KRNNAGRQQI KAAYLQETGK PLDETLKKAL TGHLEEVVLA LLKTPAQFDA DELRAAMKGL GTDEDTLIEI LASRTNKEIR DINRVYREEL  
KRD LAKDITS DTSGDFRNAL LSLAKGDRSE DFGVNEDLAD SDARALYEAG ERRKGTDVNV FNTILTTRSY PQLRRVFQKY TKYSKHD MNK  
VLDLELKGDI EKCLTAIVKC ATSKPAFFAE KLHQAMKGVG TRHKALIRIM VSRSEIDMND IKAFYQKMYG ISLCQAILDE TKG DY EKILV ALCGGN

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human Annexin A1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.