

Human SEMG1 Recombinant Protein

Catalog Number: PROTP04279

Overview

Product Name	Human SEMG1 Recombinant Protein
Description	Human SEMG1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 24-462aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 7.5) containing 0.5M NaCl, 10% glycerol, 250mM Imidazole, 0.1mM PMSF, 1mM DTT

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, P04279, 24-462aa; MGSSHHHHHH SGLVPRGSH MGSQKGGSKG RLPSEFSQFP HGQKGQHYSG QKGKQQTESK GSFSIQYTYH VDANDHDQSR KSQQYDLNAL HKTTKSQRHL GGSQQLLHMK QEGRDHDKSK GHFHRVVIHH KGGKAHRGTQ NPSQDQGNP SGKGISSQYS NTEERLWVHG LSKEQTSVSG AQKGRKQGGG QSSYVLQTEE LVANKQQRET KNSHQNKGHY QNVVEVREEH SSKVQTSLCP AHQDKLQHGS KDIFSTQDEL LVYNKNQHQT KNLNQDQHQG RKANKISYQS SSTEERLLHY GENGVQKDVS QSSIYSQTEE KAQGKSQKQI TIPSQEQEHS QKANKISYQS SSTEERLLHY GENGVQKDVS QRSIYSQTEK LVAGKSIQA PNPKQEPWHG ENAKGESGQS TNREQDLLSH EQKGRHQHGS HGGLDIVIIE QEDDSRHLA QHLNDRNPL FT

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human SEMG1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.