

Mouse IL-2 (C160S) Recombinant Protein

Catalog Number: PROTP04351

Overview

Product Name	Mouse IL-2 (C160S) Recombinant Protein
Description	Mouse IL-2 (C160S) Recombinant Protein expressed in E. coli (non-tagged). Sequence domain: 21-169aa. Application(s): SDS-PAGE, Bioactivity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	Non-Tagged
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) 10% glycerol

Concentration

0.5mg/ml (determined by BCA assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Mouse, P04351, 21-169aa; MAPTSSSTSS STAEAQQQQQ QQQQQQHLE QLLMDLQELL SRMENYRNLK LPRMLTFKfy LPKQATELKD LQCLEDELGP LRHVLDLTQS KSFQLEDAEN FISNIRVTVV KLKGSNTFE CQFDDESATV VDFLRRWIAF SQSIISTSPQ

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse IL-2 (C160S) Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.