

## Human ErbB2/Her2 Recombinant Protein

Catalog Number: PROTP04626

### Overview

Product Name	Human ErbB2/Her2 Recombinant Protein
Description	Human ErbB2/Her2 Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 23-652aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

0.25mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 85% by SDS-PAGE

### Amino Acid Sequence

Human, P04626, 23-652aa; TQVCTGDMK LRLPASPETH LDMLRHLYQG CQVVQGNLEL TYLPTNASLS FLQDIQEVQG YVLIAHNQVR QVPLQRLRIV RGTQLFEDNY ALAVLDNGDP LNNTPVTGA SPGGLRELQL RSLTEILKGG VLIQRNPQLC YQDTILWKDI FHKNNQLALT LIDTNRSRAC HPCSPMCKGS RCWGESSEDC QSLTRTV CAG GCARCKGPLP TDCCHEQCAA GCTGPKHSDC LACLHFNHSG ICELHCPALV TYNTDTFESM PNPEGRYTFG ASCVTACPYN YLSTDVGSCT LVCPLHNQEV TAEDGTQRCE KCSKPCARVC YGLGMEHLRE VRAVTSANIQ EFAGCKKIFG SLAFLPESFD GDPASNTAPL QPEQLQVFET LEEITGYLYI SAWPDSLPLD SVFQNLQVIR GRILHNGAYS LTLOQLGISW LGLRSLRELG SGLALIHHT HLFCVHTVPW DQLFRNPHQA LLHTANRPED ECVGEG LACH QLCARGHCWG PGPTQCVNCS QFLRGQECVE ECRVLOGLPR EYVNRHCLP CHPECQPQNG SVTCFGPEAD QCVACAHYKD PPFVCVARCPS GVKPDL SYMP IWKFPDEEGA CQPCPINCTH SCVDLDDKGC PAEQRASPLT

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human ErbB2/Her2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.