

Rat Lactate Dehydrogenase A/LDHA Recombinant Protein

Catalog Number: PROTP04642

Overview

Product Name	Rat Lactate Dehydrogenase A/LDHA Recombinant Protein
Description	Rat Lactate Dehydrogenase A/LDHA Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 1-332aa. Application(s): Enzyme Activity, SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 1mM DTT, 20% glycerol

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Rat, P04642, 1-332aa; MAALKDQLIV NLLKEEQVPQ NKITVVGVA VGMACAISIL MKDLADELAL VDVIEDKCLKG EMMDLQHGS LFLKTPKIVSS KDYSVTANSK LVIITAGARQ QEGESRLNLV QRNVNIFKFI IPNVVKYSPQ CKLLIVSNPV DILTYVAWKI SGFPKNRVIG SGCNLD SARF RYLMGERLGV HPLSCHGWVL GEHGDSSVPV WSGVNVAGVS LKSLNPQLGT DADKEQWKDV HKQVVD SAYE VIKLKG YTSW AIGLSVADLA ESIMKNLRRV HPITMIKGL YGIKEDVFLS VPCILGQNGI SDVVKVTLTP DEEARLK KSA DTLWGIQKEL QFLEHHHHHH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Rat Lactate Dehydrogenase A/LDHA Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.