

Human IGF-I Recombinant Protein

Catalog Number: PROTP05019

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human IGF-I Recombinant Protein
Description	Human IGF-I Recombinant Protein expressed in E. coli (non-tagged). Sequence domain: 49-118aa. Application(s): SDS-PAGE, Bioactivity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 20ug
Tag	Non-Tagged
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4)

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P05019, 49-118aa; MGPETLCGAE LVDALQFVCG DRGFYFNKPT GYGSSRRAP QTGIVDECCF RSCDLRRLEM YCAPLKPAKS A

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human IGF-I Recombinant Protein](#)

For Research Use Only. Not for use in diagnostic procedures.