

Human c-Jun Recombinant Protein

Catalog Number: PROTP05412

Overview

Product Name	Human c-Jun Recombinant Protein
Description	Human c-Jun Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-241aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 2mM DTT, 0.1M NaCl.

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P05412, 1-241aa; MTAKMETTFY DDALNASFLP SESGPYGYSN PKILKQSMTL NLADPVGSLK PHLRAKNSDL LTSPDVGLLK LASPELERLI IQSSNGHITT TPTPTQFLCP KNVTDQEGF AEGFVRALAE LHSQNTLPSV TSAAQPVNGA GMVAPAVASV AGGSGSGGFS ASLHSEPPVY ANLSNFNPGA LSSGGGAPSY GAAGLAFPAQ PQQQQPPHH LPQQMPVQHP RLQALKEEPQ TVPEMPGETP P

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human c-Jun Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.