

Mouse Renin Recombinant Protein

Catalog Number: PROTP06281

Overview

Product Name	Mouse Renin Recombinant Protein
Description	Mouse Renin Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 22-402aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 4ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Mouse, P06281, 22-402aa; LPTRTATFER IPLKKMPSVR EILEERGVDMLRLSAEWGVF TKRPSLTNLT SPVVLTYLN TQYYGEIGIG TPPQTFKVI
DTGSANLWVP STKCSRLYLA CGIHSLYESS DSSSYMENGS DFTIHYGSGR VKGFLSQDSV TVGGITVTQT FGEVTEPLI PFMLAKFDGV
LGMGFPAQAV GGVTPVFDHI LSQGVLEKEV FSVYYNRGSH LLGGEVVLGG SDPQHYQGNF HYVSISKTDS WQITMKGVSV GSSTLLCEEG
CAVVVDTGSS FISAPTSSLK LIMQALGAKE KRIEYVVCN SQVPTLPDIS FDLGGRAYTL SSTDYVLQYP NRRDKLCTLA LHAMDIPPT
GPVWVLGATF IRKFYTEFDR HNNRIGFALA R

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Renin Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.