

## Mouse Cathepsin L Recombinant Protein

Catalog Number: PROTP06797

### Overview

Product Name	Mouse Cathepsin L Recombinant Protein
Description	Mouse Cathepsin L Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 18-334aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 4ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

0.25mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Mouse, P06797, 18-334aa; TPKFDQTFSA EWHQWKSTHR RLYGTNEEEW RRAIWEKNMR MIQLHNGEYS NGQHGFSEMEM NAFGDMTNEE  
FRQVVNGYRH QKHKKGRLFQ EPLMLKIPKS VDWRKGCVT PVKNQGCQS CWAFSASGCL EGQMFLKTGK LISLSEQNLV DCSHAQGNQG  
CNGGLMDFAF QYIKENGGLD SEESYPYEAQ DGSCKYRAEF AVANDTGFVD IPQKEKALMK AVATVGPISV AMDASHPSLQ FYSSGIYYEP  
NCSSKNLDHG VLLVGYGYEG TDSNKNKYWL VKNSWGSEWG MEGYIKIAKD RDNHCGLATA ASYPVVN

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Mouse Cathepsin L Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.