

## Human Hexosaminidase A/HEXA Recombinant Protein

Catalog Number: PROTP06865

### Overview

Product Name	Human Hexosaminidase A/HEXA Recombinant Protein
Description	Human Hexosaminidase A/HEXA Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 23-529aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4)

### Concentration

0.5mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, P06865, 23-529aa; LWPWPQNFQT SDQRYVLYPN NFQFQYDVSS AAQPGCSVLD EAFQRYRDLL FGSGSWPRPY LTGKRHTLEK  
NVLVSVVTP GCNQLPTLES VENYTLTIND DQCLLSETV WGALRGLET FSQLVWKAEG TFFINKTEIE DFPRFPHRGL LLDTSRHYLP LSSILDTLDV  
MAYNKLNVFH WHLVDDPSFP YESFTFPELM RKGSYNPVTH IYTAQDVKEV IEYARLRGIR VLAEFDTPGH TLSWGPGIPG LLTPCYSGSE  
PSGTFGPVNP SLNNTYEFMS TFFLEVSSF PDFYLHLGGD EVDFTCWKSN PEIQDFMRKK GFGEDFKQLE SFYIQTLLDI VSSYGKGYVV  
WQEVFDNKVK IQPDTIIQVW REDIPVNYMK ELEVTKAGF RALLSAPWYL NRISYGPDWK DFYIVEPLAF EGTPEQKALV IGGEACMWGE  
YVDNTNLVPR LWPRAGAVAE RLWSNKLTSD LTFAYERLSH FRCELLRRGV QAQPLNVGFC EQEFEQT

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Hexosaminidase A/HEXA Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.