

KLK1 Kallikrein-1 Human Recombinant Protein

Catalog Number: PROTP06870

Overview

| | |
|--------------|---|
| Product Name | KLK1 Kallikrein-1 Human Recombinant Protein |
| Description | KLK1 Human Recombinant produced in HEK293 cells is a single, glycosylated, polypeptide chain containing 244 amino acids (Pro19-Ser262). KLK-1 is expressed with an 6 amino acid His tag at C-Terminus and purified by proprietary chromatographic techniques. |
| Size | 2ug, 10ug, 1mg |
| Tag | |
| Form | Sterile Filtered clear solution. |
| Source | HEK293 |
| Formulation | Provided as a solution in 20mM TrisHCl, 150mM NaCl, 2mM CaCl2, pH- 8. |

Concentration

Storage

Store KLK1 at 4°C if entire vial will be used within 2-4 weeks.Store frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please avoid freeze thaw cycles.

Purity

Greater than 95.0% as determined by SDS-PAGE.

!description_after_attributes! ?/description_after_attributes?

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



KLK1 Kallikrein-1 Human Recombinant Protein