

Human Glutathione Peroxidase 1/GPX1(U49C) Recombinant Protein

Catalog Number: PROTP07203

Overview

Product Name	Human Glutathione Peroxidase 1/GPX1(U49C) Recombinant Protein
Description	Human Glutathione Peroxidase 1/GPX1(U49C) Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-203aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 2mM DTT, 30% glycerol, 100mM NaCl

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P07203, 1-203aa; MCAARLAAAA AAAQSVYAFS ARPLAGGEPV SLGSLRGKVL LIENVASLCG TTVRDYTMN ELQRRLGPRG
LVVLGFPCNQ FGHQENAKNE EILNSLKYVR PGGGFEPNFM LFEKCEVNGA GAHPLFAFLR EALPAPSSDA TALMTDPKLI TWSPVCRNDV
AWNFEKFLVG PDGVPLRRYS RRFQTIDIEP DIEALLSQGP SCA

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Glutathione Peroxidase 1/GPX1(U49C) Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.