

Human recombinant Thrombomodulin/BDCA-3 protein

Catalog Number: PROTP07204-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human recombinant Thrombomodulin/BDCA-3 protein
Description	Human recombinant Thrombomodulin/BDCA-3 protein expressed in Baculovirus with His-Tag. Sequence domain: 22-515aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P07204, 22-515aa; EPQPGGSQCV EHD CFALYPG PATFLNASQI CDGLRGHLMT VRSSVAADVI SLLNGDGGV GRRRLWIGLQ LPPGCGDPKR LGPLRGFQWV TGDNNTSYSR WARLDLNGAP LCGPLCVAVS AA EATVPSEP IWEEQQCEVK ADGFLCFEHF PATCRPLAVE PGAAAAVSI TYGTPFAARG ADFQALPVGS SAAVAPLGLQ LMCTAPPGAV QGHWAREAPG AWDCSVENGG CEHACNAIPG APRCQCPAGA ALQADGRSCT ASATQSCNDL CEHFCVPNPD QPGSYSCMCE TGYRLAADQH RCEDVDDCIL EPSPCPQRCV NTQGGFECHC YPNYDLVDGE

CVDPVPCFRANCEYQCQPLNQTSYLCVCAEGFAPIPEPHRCQMFCNQTACPADCDPNTQASCECPEGYILDDGFICTDIDECENGGFC
SGVCHNLPGTFECICGPDSA LARHIGTD CD SGKVDGGDSGSGEPSPSTPGSTLTTPAVGLVHS

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human recombinant Thrombomodulin/BDCA-3 protein](#)

For Research Use Only. Not for use in diagnostic procedures.