

Mouse recombinant SPARC protein

Catalog Number: PROTP07214

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Mouse recombinant SPARC protein
Description	Mouse recombinant SPARC protein expressed in Baculovirus with His-Tag. Sequence domain: 18-302aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by Absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Mouse, P07214, 18-302aa; APQQTEVAEE IVEEETVVEE TGVPVGANPV QVEMGEFEDG AEETVEEVVA DNPCQNHHC HGKVCELDES
NTPMCVCQDP TSCPAPIGEF EKVCSNDNKT FDSSCHFFAT KCTLEGTKKG HKLHLDYIGP CKYIAPCLDS ELTEFPLMR DWLKNVLVTL
YERDEGNNLL TEKQKLRVKK IHENEKRLA GDHPVELLAR DFEKNYNMYI FPDHWQFGQL DQHPIDGYLS HTELAPLRAP LIPMEHCTTR
FFETCDLDND KYIALEEWAG CFGIKEQDIN KDLVI

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Mouse recombinant SPARC protein](#)

For Research Use Only. Not for use in diagnostic procedures.