

Human recombinant PROS1 protein, AF

Catalog Number: PROTP07225

Overview

Product Name	Human recombinant PROS1 protein, AF
Description	PROS1, also known as protein S, is a vitamin K-dependent protein that plays a crucial role in the regulation of blood coagulation. It functions as a cofactor for protein C, another key player in the coagulation cascade. PROS1 acts as an anticoagulant by inhibiting the activity of certain clotting factors and promoting the breakdown of blood clots. It helps maintain the balance between clotting and anticoagulation, preventing the formation of excessive blood clots that can lead to thrombosis.
Size	5ug,20ug,100ug
Tag	His-SUMO Tag (N-term)
Form	Lyophilized
Source	Escherichia coli
Formulation	The protein was lyophilized from a 0.2 um filtered solution containing 1X PBS, pH 7.4. If you have any concerns or special requirements, please confirm with us.

Storage

Lyophilized protein should be stored at -20°C for 1 year. Upon reconstitution, store at 2°C to 8°C for up to 1 week. Further dilute in a buffer containing a carrier protein or stabilizer (e.g. 0.1% BSA, 10%FBS, 5%HSA or 5% trehalose solution), protein aliquots should be stored at -20°C or -80°C for 3-6 months. Avoid repeated freeze/thaw cycles.

Purity

>95% as determined by SDS-PAGE.

Molecular weight

The protein has a calculated MW of 82.89 kDa. The protein migrates as 75 kDa under reducing condition (SDS-PAGE analysis).

Activity

Testing in process

Endotoxin

<1 EU per 1 ug of the protein by the LAL method.

Amino Acid Sequence

A DNA sequence encoding Human PROS1 Protein (#P07225)(Ala42-Ser676) was expressed with polyhistidine-SUMO tag at the N-terminus.

Reconstitution

Centrifuge at 3000 rpm for 5 mins before opening. It is recommended to reconstitute the lyophilized protein in sterile H₂O to a concentration not less than 100 ug/mL and incubate the stock solution at room temperature for at least 20 mins to ensure sufficient re-dissolved. Do Not Vortex!

Vigorous shaking may impair the biological activity of the protein.

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human recombinant PROS1 protein, AF

For Research Use Only. Not for use in diagnostic procedures.