

Mouse recombinant Erythropoietin/EPO Protein

Catalog Number: PROTP07321

Overview

Product Name	Mouse recombinant Erythropoietin/EPO Protein
Description	Mouse recombinant Erythropoietin/EPO Protein protein expressed in Baculovirus with His-Tag. Sequence domain: 27-192aa. Application(s): SDS-PAGE, Bioactivity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by Absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Mouse, P07321, 27-192aa; MAPPRLI CDSRVLERYI LEAKEAENVT MGCAEGPRLS ENITVPDTKV NFYAWKRMEV EEQAIEVWQG LLLSEAILQ
AQALLANSSQ PPETLQLHID KAISGLRSLT SLLRVLGAQK ELMSPDTP PAPLRTLVD TFCKLFRVYA NFLRGKCLKY TGEVCRRGDR

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse recombinant Erythropoietin/EPO Protein

For Research Use Only. Not for use in diagnostic procedures.