

Human M-CSF R/CD115 Recombinant Protein

Catalog Number: PROTP07333

Overview

Product Name	Human M-CSF R/CD115 Recombinant Protein
Description	Human M-CSF R/CD115 Recombinant Protein expressed in Baculovirus with hIgG-His-tag. Sequence domain: 20-517aa. Application(s): SDS-PAGE, Bioactivity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hIgG-His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P07333, 20-517aa; IPVIEPSVPE LVVKPGATVT LRCVGNNGSVE WDGPPSPHWT LYSDGSSSIL STNNATFQNT GTYRCTEPGD
PLGGSAAIHL YVKDPPARPWN VLAQEVVVE DQDALLPCLL TDPVLEAGVS LVRVRGRPLM RHTNYSFSPW HGFTIHRAKF IQSQDYQCSA
LMGGRKVMMSI SIRLKVQKVI PGPPALTLVP AELVRIRGEA AQIVCSASSV DVNFDVFLQH NNTKLAIPQQ SDFHNNRYQK VLTLNLDQVD
FQHAGNYSVCV ASNVQKGHST SMFFRVVESA YLNLSSSEQNL IQEVTVGEGL NLKVMVEAYP GLQGFNWTYL GPFSDHQPEP KLANATTKDT
YRHTFTLSLP RLPKPEAGRY SFLARNPGGW RALTFELTLR YPPEVSVIWT FINGSGTLLC AASGYQPNV TWLQCSGHTD RCDEAQLVQV
WDDPYPEVLS QEPFHKVTVQ SLLTVETLEH NQTYECRAHN SVGSGSWAFI PISAGAHTHP PDEFLFTP

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human M-CSF R/CD115 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.