

Human HSP90 alpha Recombinant Protein

Catalog Number: PROTP07900

Overview

Product Name	Human HSP90 alpha Recombinant Protein
Description	Human HSP90 alpha Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-732aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Liquid. In Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P07900, 1-732aa; MPEETQTQDQ PEEEEVETAFQAEIAQLMSLIINTFYNSKEIFLRELISNSSDALDKIRYESLTDPSKLDGKELHINLIPNKQDRTLTIVDTGIGMTKADLNNLGTIAKSGTKAFMEALQAGADISMIGQFGVGFYSAYLVAEKVTVITKHNDDEQYAWESSAGGSFTVRTDTGPEMGRGTKVILHLKEDQTEYLEERRIKEIVKKHSQFIGYPITLVFEKERDKEVSDDEAEEKEDKEEKEKEEKESEDKPEIEDVGSDEEEEEKDGDGKKEKKEKKEKYIDQEELNKTPIWTRNPDDITNEEYGEFYKSLTNDWEDHLAVKHFSVEGQLEFRALFVPRRAPDFLFENRKKKNNIKLYVRRVFI MDNCEELIPEYLNFI RGVVD SEDLPLNISREMLQQSKILK VIRKNLVKKLLELFTELAEDKENYKFFYEQFSKNIKLGIIHEDSQNRKKLS ELLRYYTSASGDENVSLKDYCTRMKENQKH IYYITGETKDVANSFAFVERLRKHGLEVIYMIPEIDCYCVQLKFEFGKTLVSVTKEGLELPEDEEEKKEEKEKTKFENLCKIMKDILEKKVEKVVVSNRLVTSPCCIVTSTYGWTANMERIMKAQALRDNSTMGYMAAKKHLEINPDHSIIETLRQKAEADKNDKSVKDLVILLYETALLSSGFSLEDPQTHANRIYRMIKLG LGIDE DDPTADDTSA AVTEEMPPLE GDDDTSRMEEVD

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human HSP90 alpha Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.