

## Human HNRNPC Recombinant Protein

Catalog Number: PROTP07910

### Overview

Product Name	Human HNRNPC Recombinant Protein
Description	Human HNRNPC Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-293aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 5mM DTT, 30% glycerol

### Concentration

0.25mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 85% by SDS-PAGE

### Amino Acid Sequence

Human, P07910, 1-293aa; MGSSHHHHHH SGLVPRGSH MASNVTNKTD PRSMNSRVFI GNLNLTLVVKK SDVEAIFSKY GKIVGCSVHK  
GFAFVQYVNE RNARA AVAGE DGRMIAGQVL DINLAAEPKV NRGKAGVKRS AAEMYGSSFD LDYDFQRDYY DRMYSYPARV PPPPIARAV  
VPSKRQRVSG NTSRRGKSGF NSKSGQRGSS KSGKLGDDL QAIKKELTQI KQKVDLLEN LEKIEKEQSK QAVEMKNDKS EEEQSSSVK  
KDETNVKMES EGGADDSAE E GDLLDDDDNE DRGDDQLELI KDDEKEAEEG EDDRDSANGE DDS

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human HNRNPC Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.