

Rat IGF-1 Recombinant Protein

Catalog Number: PROTP08025

Overview

Product Name	Rat IGF-1 Recombinant Protein
Description	Rat IGF-1 Recombinant Protein expressed in E. coli (non-tagged). Sequence domain: 49-118aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	Non-Tagged
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Rat, P08025, 49-118aa; MGPETLCGAE LVDALQFVCG PRGFYFNKPT GYGSSIRRAP QTGIVDECCF RSCDLRRLEM YCAPLKPTKS A

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Rat IGF-1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.